

# **BI+S1**

**Digital Endoscope**



[www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)

**Quality, Professional, Innovative**

# Handheld Digital Endoscope S1

## Specifications

### Control Unit

Resolution	720 x 480 pixels
Light source	Natural white (LED)
LCD screen	3.5 inch colorful TFT-LCD
Magnification	8~20 times
Image format	JPEG (photo) , H.264 (video)
Interface	Mini USB port, AV out port
Media	Media Micro SD flash card, 1G
Power supply	3.7V / 2200mAh rechargeable lithium battery
External power	AC 100~240 V, 50/60 Hz
Working time	3 hours
Charging time	5 hours
Package	280mm (L)*220mm (W) *140mm(H)

### Our Images



### Lens

Digital Otoscope (ear inspection)	
Focus range	10~50 mm
Diameter	5.5mm
Length	70mm
Weight	31g



Digital Laryngoscope (throat inspection)	
Focus range	80~160mm
Diameter	230mm
Length	70mm
Weight	62g



Digital Rhinoscope (nasal inspection)	
Focus range	10~50 mm
Diameter	5.5mm
Length	250mm
Weight	31g



**Wuxi Biomedical Technology Co., Ltd.**

[www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: [sales@cnbiomedical.com](mailto:sales@cnbiomedical.com)

Mobile phone: 86-13382885910

**BIO** is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.